Elite ORP Tester Operating Instructions

68X001318 Rev 0 June, 2017

English

This instruction manual is applicable for Elite ORP Tester.

Installation and Uninstallation of the Battery

The Elite ORP tester uses four AAA 1.5 V batteries.

1. To remove the battery cover, See Figure 1. Clear the front catch and then the back catch, before sliding the cover off.
2. To remove the battery plate, push the center tab towards the front of the tester as show in Figure 2. Once unlocked, remove the plate to access the batteries.
3. Invert the tester upside down to remove the batteries. Each side uses two AAA batteries. Orient each battery with positive terminal facing downward.
4. To lock the battery place, align the small tabs into the guide ribs on the housing and then press down. See Figure 4.


Figure 1: Removing battery cover


Figure 5: Battery plate unlock position

Figure 6: Battery plate lock position

Figure 2: Push to unlock


Figure 3: Align tabs


Figure 4: Push down to lock

Getting Started

The Elite ORP tester has been factory calibrated and usually works well out of the box. However, after extended periods of non-use, it is best to remove the sensor cap and soak the sensor in warm tap water or ORP buffer for 10 minutes or so. A brief rinse with deionized (DI) water is OK, but avoid soaking or storing in deionized water as this will shorten ORP electrode life. Prior to taking the measurements, periodic calibration with certified standards is recommended for best accuracy.


Temperature Settings

1. Press to enter setup window. Press to select Temp set. The display shows Set °C/°F and Temp Cal.
2. Press to select Set °C/°F. The display automatically shows °C and °F.
3. Scroll down by pressing to toggle between °C and °F.
4. Press to select a temperature unit. The display shows the selected temperature setting with a ✓.

Backlight Settings


1. Press to enter setup window. Scroll down by pressing to select Settings. Press to select Settings. The display shows Backlight.
2. Press to select Backlight. The display shows ON and OFF.
3. Scroll down by pressing to toggle between ON and OFF. Backlight ON increase readability in low light conditions.
4. Press to select the desired backlight option. The display shows the selected backlight option with a .

ORP Calibration

Calibration should be done regularly. Preferred once in a week. ORP testers allow one point calibration.

1. Press to power on if needed.
2. Dip electrode about 2 cm to 3 cm in the calibration standard solution.
3. Press the to enter calibration mode.
4. Allow the reading to stabilize. The timer icon flashes during the reading stabilization. Once the reading stabilized the timer stops flashing.
5. Upper display shows current measured value based on the factory default.
6. Lower display shows current measured mV with last calibration (if no calibration performed earlier both display shows the same measurement value).
7. Press to adjust the mV value to match the calibration standard value.
8. Press to confirm the calibration or press key to exit the calibration without confirmation.

ORP Measurement

1. Press to power on if needed.
2. Dip the electrode about 2 cm to 3 cm into the test solution. Stir and let the reading stabilize. The timer icon flashes during this time. Once the reading is stabilized, the timer stops flashing and will appear to indicate the stability of the reading.

CAUTION: Testing dry samples is not accurate and can lead to sensor damage or breakage. Soils must be wet and free of particulates that may scratch the glass sensor. Excessive force into dry samples can cause glass breakage.

3. Note the ORP value or press to freeze the reading. To release the reading, press again.
4. Press to turn off tester. If the key is not pressed for 8.5 minutes, the tester will automatically shut off to conserve batteries.

User Reset

Reset the ORP calibration to the default using the user reset function.

1. Press to enter setup window. Scroll down by pressing to select Reset. Press to select reset. The display shows User Reset and Factory Reset.
2. Press to select User Reset. The display automatically shows No and Yes.
3. Scroll down by pressing to toggle between No and Yes.
4. Press to confirm either No or Yes. The display shows the user reset option with a .

Factory Reset

Reset the tester to the factory default using the factory reset function.

1. Press to enter setup window. Scroll down by pressing to select Reset. Press to select reset. The display shows User Reset and Factory Reset.
2. Scroll down by pressing to toggle between the resets. Press to select Fact Reset. The display automatically shows No and yes.
3. Scroll down by pressing to toggle between No and Yes.
4. Press to confirm either No or Yes. The display shows the factory reset option with a .


Manual Temperature Setting

The MTS range is 0°C to 60.0°C (32.0°F to 140.0°F). User reset will set temperature to default value 25°C or 77°F.

1. Press to select setup window. Scroll down by pressing to select Temp Set. Press to select Temp set. The display shows Set °C/F and Temp Cal.
2. Scroll down by pressing to toggle between Set °C/F and Temp Cal. Press to select Temp Cal.
3. The lower display shows the current measured temperature reading based on the last set offset and the upper display shows the current measured temperature reading based on factory default calibration.
4. Dip the tester into a solution of known temperature and allow time for the in built temperature sensor to stabilize.
5. Press to adjust the temperature value or press the to confirm the calibrated value as new temperature value of the solution.


Note: To exit without confirming the calibration, press .


HOLD Function

This feature helps to freeze the display for a delayed observation.

1. Press  button to freeze the measurement.

2. Press  again to release the measurement.

Sensor Maintenance


1. Rinse the sensor with sensor storage solution after each measurement. Care has to be taken not to damage the sensor's glass electrode.
2. In aggressive chemicals, dirty or viscous solutions, and solutions with heavy metals or proteins, take readings quickly and rinse electrode immediately afterward.
3. If possible, keep a small piece of paper or sponge in the sensor cap – moistened with clean water or electrode storage solution (NOT de-ionized water) – and close the cap over the sensor.

Sensor Replacement

Replace the sensor module at the fraction of the cost of a new tester. When the tester fails to calibrate or gives fluctuating readings in calibration standards, need to change the electrode.

1. With dry hands, grip the ring with sensor facing you. Twist the ring clockwise. Save the ring for later use.
2. Pull the old sensor module away from the tester.
3. Align the four tabs on the new module so that they match the four slots on the tester.
4. Gently push the module onto the slots to sit it in position. Push the smaller O-ring fully onto the new sensor module. Push the other O-ring over the module and thread it into place by firmly twisting counter clockwise.

Note: It is necessary that recalibrate the tester prior to measurement after a sensor replacement.

Applications

The Elite ORP tester can be used for the following:

- Aquaculture
- Food sectors and more
- Steam generators
- Aquariums Studies
- Hydroponics
- Studies
- Boilers
- Labs
- Water and wastewater treatment
- Car washes
- Sanitation plants
- Water quality testing in pools
- Ecology
- Spas

Warranty

This instrument is supplied with a warranty against manufacturing defects for a period of one year from the date of purchase.

Return of Items

Authorization must be obtained from distributor before returning items for any reason. When applying for authorization, please include information regarding the reason the item(s) are to be returned.

We reserve the right to make improvements in design, construction and appearance of products without notice. Prices are subject to change without notice.

Specification	Elite ORP
ORP/RMV	
ORP Range	-1000 to 1000 mV
Resolution	+/- 1 mV
Relative accuracy	+/-2mV +/- 1LSB
Calibration points	One
Calibration window	+/- 200mV
Calibration type	Off Set
Temperature	
Temperature range	0°C to 60°C 32.0°F to 140.0°F
Temperature resolution	0.1°C / 0.1°F
Temperature accuracy	NA
Temperature compensation	Manual
Temperature calibration window	"0°C to 60.0°C 32.0°F to 140.0°F"

Specification	Elite ORP
General Specification	
Display	Graphics, Dot Matrix 80X100 Pixel
Backlight	Yes, Selectable (30sec from the last key press)
Auto off	8.5 minutes (From Last Key press)
Reset	User / Factory
Hold Function	Yes
Power requirement	Four AAA 1.5 V batteries
Battley life	>250 hours
Water proofing	IP67
Regulatory Certifications	CE, FCC
Environmental Operating Conditions	
Ambient Operating Temperature	5°C to 45°C
Relative Humidity	5% to 85% Non-Condensing
Storage Temperature	-20°C to 60°C
Storage Humidity	5% to 85%, Non-Condensing

Accessories

Ordering Code	Product Description
ORP Tester	
ELITEORP	Elite ORP tester with batteries
Storage and Cleaning Solutions	
ECRE005	Storage solution for ORP sensor, 480 mL
ECDPCBT	Protein remover solution, 480 mL
Other Accessories	
ELITECAP	Replacement sensor cap
TSLANYARD	Tester lanyard

Self Diagnostic Messages

	Batteries are weak and need replacement soon.
stable error	Appears when calibration is attempted but the reading is not yet stable. Wait for the reading to stabilize or manually confirm the calibration by pressing enter.
over range	The reading is above the measuring range of tester.
under range	The reading is below the measuring range of tester.

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