UNPACKING

Carefully remove your Scout balance and each of its components from the package. Store the packaging material to ensure safe storage and future transport.



Included Components:





- Balance
- Quick Start Guide
- Warranty Card
- AC Adapter & Plug
- Pan
- Calibration Mass (only for SJX/E models with 64g capacity)
- Scoop (only for models with 64g capacity)
- Capacity label set (SJX323N/E and SJX1502N/E only)

SELECTING THE LOCATION

Place on a firm, steady surface with sufficient space. Avoid excessive air currents, vibrations, heat sources, or rapid temperature changes.

PREPARATION FOR USE

Step-1	Release the red Transportation Lock on subplatform of the balance	
	<p>Non SJXN/E models: Turn the pointer 90° counter-clockwise.</p>	 <p>Locked Position</p> <p>Unlocked Position</p>
	<p>SJXN/E models: Turn the pointer 90° clockwise, then remove the shipping lock by pulling it out before installing the pan.</p> <p>Store the shipping lock properly for future transport. To reinstall it, follow the removal process in reverse order.</p> <p>Pointer pulled out:</p> 	 <p>Locked Position</p> <p>Unlocked Position</p>

<p>Step-2</p>	<p>Install the Pan</p> <p>For SJX/E: <i>Rectangular pan:</i> Place diagonally on sub-platform and rotate counter-clockwise until it locks.</p> <p><i>Round pan:</i> Set straight down on subplatform.</p> <p>For SJX: Install the stainless steel pan directly on the subplatform.</p>	 <p style="text-align: center;">For SJX/E</p>
<p>Step-3</p>	<p>Attach Capacity Label</p> <p>For SJX323N/E and SJX1502N/E: Attach the capacity label above the screen at the position as demonstrated in the graphic, and align the edges with the dotted lines.</p>	 <p style="text-align: center;">For SJX323N/E and SJX1502N/E</p>
<p>Step-4</p>	<p>Connect Power Connect the AC Adapter to the balance input jack.</p> <p><i>Battery Installation</i> Install the Four “AA” (LR6) batteries (not included) with polarity as shown in the battery compartment.</p>	 <p style="text-align: center;">AC Adapter Jack Battery Compartment</p>
<p>Step-5</p>	<p>Level the balance Level the balance before using by adjusting the feet until the air bubble is centered</p>	 <p style="text-align: center;">YES NO</p>

SAFETY PRECAUTIONS

CAUTION: Read all safety warnings before installing, making connections, or servicing this equipment. Failure to comply with these warnings could result in personal injury and/or property damage. Retain all instructions for future reference.


- Verify that the local AC power supply voltage is within the input voltage range printed on the AC adapter's ratings label.
- Only connect the AC adapter to a compatible grounded socket.
- Position the instrument such that the AC adapter can be easily disconnected from the socket.
- Position the power cord so that it does not pose a potential obstacle or tripping hazard.
- Operate the equipment only under ambient conditions specified in the user instructions.
- This equipment is intended for indoor use and should only be operated in dry locations.
- Do not operate the equipment in hazardous or explosive environments.
- Disconnect the equipment from mains power before cleaning or servicing.
- Service should only be performed by authorized personnel.


Intended Use

This instrument is intended for use in laboratories, schools, businesses and light industry. It must only be used for measuring the parameters described in these operating instructions. Any other type of use and operation beyond the limits of technical specifications, without written consent from OHAUS, is considered as not intended. This instrument complies with current industry standards and the recognized safety regulations; however, it can constitute a hazard in use. If the instrument is not used according to these operating instructions, the intended protection provided by the instrument may be impaired.

CONTROLS

- **Button Functions**

Button	 Yes	 No	 Back	 Exit
Short Press	Zero/On Turns the balance on; If balance is On, sets Zero	Print Initiates print	Function Initiates an application mode	Tare Tare
Long Press Extended Press	Off Turns the balance Off None	Units Changes the weighing unit None	Mode Changes the application mode None	Menu Enters the user menu View the preset Tare value
Short Press (Menu Function)	Yes Accepts the current setting	No Advances to the next menu or setting	Back Moves to previous menu item	Exit Exits the user menu

CALIBRATION

When the balance is operated for the first time, a calibration is recommended to ensure accurate weighing results.

- **External Span Calibration**

With balance turned on:

Step-1	Press and hold Menu until [MENU] (Menu) is displayed.
Step-2	Release the button, the display will show [C.R.L.].
Step-3	Press Yes to accept. [SPAN] will then be shown.
Step-4	Press Yes to begin the span calibration. [--C--] will blink and the Zero reading is stored automatically.
Step-5	The specified mass blinks on the display. Place the specified span calibration mass on the pan or press No to select an alternate weight. [--C--] will blink while the reading is stored. The display will show [done] if the calibration was successful. The balance returns to the previous application mode and is ready for use.

- **Internal Calibration (not available for SJX/E models)**

With balance turned on:

Step-1	Press and hold Menu until [MENU] (Menu) is displayed.	
Step-2	Release the button, the display will show [C.R.L.].	
Step-3	Press Yes to accept. [In.CAL] will then be shown.	
Step-4	Press Yes to begin the Internal calibration. Zero reading is stored automatically.	
Step-5	Pull the InCal handle while the display show [PULL].	
Step-6	Push the InCal handle while the display show [PUSH].	
Step-7	Repeat the step-5, and step-6 one more time.	
Step-8	The display will show [done] if the calibration was successful. The balance returns to the previous application mode and is ready for use.	

OPERATION & APPLICATIONS

The balance can be configured to operate in one or more Application modes. Press and hold **Mode** key to select your Application mode.

The Scout SJX balance incorporates the following Applications:

[WEIGH]	[Count]	[PERcnt]	[CHECK]	[total]	[SPEC.Gr]	[HoLd]
Weighing	Parts Counting	Percent Weighing	Check Weighing	Totalization	Density Determination	Display Hold

- **Weighing**

Step-1	Press the Zero/On button to turn the balance on.
Step-2	Select Weighing (default application) by pressing & holding Mode key.
Step-3	Press Tare or Zero if necessary to begin.
Step-4	Place objects on the pan to display the weight. When stable, the * appears on the display. The display then shows the resulting weight of the sample.

MENU STRUCTURE

<u>Sub-Menus</u>	<u>C.a.l**</u>	<u>S.E.t.U.P</u>	<u>M.o.d.E</u>	<u>U.n.i.t**</u>	<u>E.n.d</u>
User Menu Items:	In.CAL *	Reset	Reset	Reset	
 Long press – Enter menu	CAL.Adj* -100	Yes/No	Yes/No	g	
 Short press - Accept	.. 0	Auxiliary	Weighing	kg	
 Short press – Go Next/Prev. menu	.. 100	Graduation**	On/Off	ct	
 Short press – Exit menu	Span	On/Off	Counting	oz	
	Lin	Filter	On/Off	ozt	
		Low/Med/High	Percent	lb	
		AZT	On/Off	gm	
		0.5/1/3/Off	Check	hkt	
		StableRange	On/Off	sgt	
		0.5/1/2/5	Totalize	twt	
		Backlight	On/Off	t	
		Off/On/Auto	Density	tola/tical	
		Auto Tare	On/Off		
		Off/On/On-acc	Hold		
		Auto Off	Disp/Peak/Off		
		Off/1/5/10			
	End	End	End	End	

Note: Bold fonts are default settings in the menus.

* These sub-menus are not available in SJX/E models.

** Available Settings may vary by local regulations.


ADDITIONAL FEATURES

- **Connecting the Interface**

Use an optional interface connectivity kit to connect the balance either to a computer, printer or OHAUS auxiliary display. Interface connection on the rear of the balance:

Below Interface kit accessories are available:


RS232, USB Host, USB Device, Ethernet, Bluetooth.

* Interface kits may vary according to local regulations

- **Weigh Below Hook**

The weigh below hook on SJX/E models is located at the reverse side of the battery cover as shown below (on SJX models it is located in the weigh below protective cover). To use this feature, remove the red protective cover underneath for the weigh below opening.

MAINTENANCE

Cleaning

WARNING: Electric Shock Hazard. Disconnect the equipment from the power supply before cleaning.

The housing may be cleaned with a cloth dampened with a mild detergent if necessary.

Attention: Do not use solvents, chemicals, alcohol, ammonia or abrasives to clean the housing or control panel.

Troubleshooting

The following table lists common problems and possible causes and remedies. If the problem persists, contact OHAUS or your authorized dealer.

Symptom	Possible Cause
Cannot turn on	No power to balance
Poor accuracy	Improper calibration Unstable environment
Cannot calibrate	Security switch in locked position
Err B.1	Weight exceeds power on zero range.
Err B.2	Weight below power on zero range.
Err B.3	Over load (weight exceeds rated capacity)
Err B.4	Under load (pan removed)
Err B.5	Tare out of range
Err B.6	Displayed value exceeds 999999 (possible in Totalization mode)
Err B.7	Unknown internal weight position
rEF.Err	Parts counting or percentage error – sample weight <1d. Balance shows error then exits parts counting or goes to [CLr.APW] .
Lo.rEF	Percent Reference weight or APW is too low for accurate results
CL E	Fail to do calibration.
USb.Err	Cannot find menu or app file in U-disk.

Service Information

If the troubleshooting section does not resolve or describe your problem, contact your authorized OHAUS service agent. For service assistance or technical support in the United States call toll-free 1-800-672-7722 ext. 7852 between 8:00 AM and 5:00 PM EST. An OHAUS product service specialist will be available to provide assistance. Outside the USA, please visit our web site, www.ohaus.com to locate the OHAUS office nearest you.

Accessories

DESCRIPTION	Item Number
RS232 kit	30268982
USB Host kit	30268983
USB Device Kit	30268984
Bluetooth Kit*	30268985
Ethernet Kit	30268986
Stacking Kit, x6	30268987

DESCRIPTION	Item Number
Stacking Kit, x1	30268988
Specific Gravity kit	30269020
Auxiliary Display Kit	30269019
Carrying Case	30269021
In-Use Cover	30269022
Printers and Cables	Contact

Note: * Bluetooth kit is only available in certain regions according to the local regulations.

LEGAL FOR TRADE (SJX...N/E models only)

When the balance is used in trade or a legally controlled application it must be set up, verified and sealed in accordance with local weights and measures regulations. It is the responsibility of the purchaser to ensure that all pertinent legal requirements are met. As the requirements vary by jurisdiction, the purchaser is advised to contact their local weights and measures office for instructions about putting the balance into service.

Settings

Before verification and sealing, perform the following steps:

1. Verify that the menu settings meet the local weights and measures regulations.
2. Units menu should be reviewed. Verify the units turned **On** are permitted by the local weights and measures regulations.
3. Perform a calibration.
4. Set the position of the Security Switch to the locked position.

Verification


A weights and measures official must perform the verification procedure. Contact the local weights and measures office for more information.

Sealing


After the Balance has been verified, it must be sealed by the weights and measures official to prevent undetected access to the legally controlled settings. Before sealing the device, ensure that the security switch is in the Locked position.

If using a wire seal, pass the sealing wire through the holes in the security switch and Bottom Housing as shown.


If using a paper seal, place the seal over the security switch and Bottom Housing as shown.

Unlocked

Locked with Paper Seal

Locked with Wire Seal

TECHNICAL DATA

Equipment Ratings:

Altitude: Up to 2000 m

Indoor use only

Operating temperature range: +5 °C to +40 °C

Specified temperature range: Refer to tables below

Relative humidity: 10% to 80% at 31°C, decreasing linearly to 50% at 40°C, non-condensing

Power: AC power adaptor input 100-240V 50/60 Hz and output 5 V DC 1 A (For use with certified or approved power supply, which must have a SELV and limited energy circuit output.), or 4 AA batteries

Main supply voltage fluctuations: up to $\pm 10\%$ of the nominal voltage

Pollution degree: 2

Installation category: II

General Models

Model	SJX323	/	SJX622	SJX1502
	SJX323/E	SJX322/E	SJX622/E	SJX1502/E
Capacity	320 ct (64 g)	320 g	620 g	1500 g
Readability	0.005 ct (0.001 g)	0.01 g	0.01 g	0.01 g
Repeatability (Std. Dev.)	0.002 g	0.01 g	0.01 g	0.02 g
Linearity	0.003 g	0.01 g	0.02 g	0.03 g
Span Calibration Mass*	50 g	200 g	300 g	1.5 kg
Linearity Calibration Mass	30, 60 g	200, 300 g	300, 600 g	1 kg, 1.5 kg
Capacity x Readability (kg)	/	/	/	1.5 x 0.00001
Capacity x Readability (ct)	320 x 0.005	1600 x 0.05	3100 x 0.05	7500 x 0.05
Capacity x Readability (oz)	2.25750 x 0.00005	11.2880 x 0.0005	21.8700 x 0.0005	52.9110 x 0.0005
Capacity x Readability (ozt)	2.05765 x 0.00005	10.2880 x 0.0005	19.9335 x 0.0005	48.2260 x 0.0005
Capacity x Readability (dwt)	41.153 x 0.001	205.76 x 0.01	398.67 x 0.01	964.52 x 0.01
Capacity x Readability (lb)	/	/	1.36690 x 0.00005	3.30690 x 0.00005
Capacity x Readability (grn)	987.66 x 0.02	4938.4 x 0.2	9568.0 x 0.2	23148.4 x 0.2
Stabilization Time (s)	1.5	1		1.5
Construction	ABS plastic housing with 304 stainless steel (SST) pan			
Draftshield	Yes		No	
Calibration	SJX/E: User-selectable external span or linearity calibration SJX: Manual Internal Calibration / User-selectable external span or linearity calibration			
Tare Range	Full Capacity by subtraction			
Weighing Units**	g, kg, ct, oz, ozt, dwt, lb, grn, hkt, sgt, twt, tical, tola			
Application Modes**	Weighing, Parts Counting, Percent Weighing, Check Weighing, Totalization, Density Determination, Display Hold			
Typical Battery Life	80 hours			
Specified Temperature	10°C (50°F) to 40°C (104°F)			
Storage Conditions	-20°C (-4°F) to 55°C (131°F) at 10% to 90% relative humidity, non-condensing			
Communication	RS232, USB Host, USB Device, Ethernet or Bluetooth*** (all available as accessory)			
Display Type	Backlit LCD: 6-digit 7-segment with white LED backlight			
Display Size	20 mm digits / 0.78 in			
Pan Size (W x D)	Ø93 mm / 3.7 in	Ø120 mm / 4.7 in		170 x 140 mm / 6.7 x 5.5 in
Scale Dimensions - SJX/E (W x D x H)	202 x 222 x 103 mm / 8 x 8.7 x 4.1 in	202 x 224 x 54 mm / 8 x 8.8 x 2.1 in		
Scale Dimensions - SJX (W x D x H)	202 x 230 x 114 mm	202 x 230 x 68 mm		
Shipping Dimensions - SJX/E (W x D x H)	300 x 250 x 129mm / 11.8 x 9.8 x 5.1 in	300 x 250 x 86 mm / 11.8 x 9.8 x 3.4 in		
Shipping Dimensions - SJX (W x D x H)	300 x 250 x 129mm			
Net Weight (SJX/E)	1 kg / 2.2 lb			
Net Weight (SJX)	1.1 kg			1.6 kg
Shipping Weight (SJX/E)	1.5 kg / 3.3 lb			
Shipping Weight (SJX)	1.7 kg			2.2 kg

Model	/	/	SJX6201	/
	SJX621/E	SJX3201/E	SJX6201/E	SJX8200/E
Capacity	620 g	3200 g	6200 g	8200 g
Readability	0.1 g	0.1 g	0.1 g	1 g
Repeatability (Std. Dev.)	0.1 g	0.1 g	0.1 g	1 g
Linearity	0.1 g	0.1 g	0.2 g	1 g
Span Calibration Mass*	300 g	3 kg	5 kg	8 kg
Linearity Calibration Mass	300, 600 g	1.5 kg, 3 kg	3 kg, 6 kg	4 kg, 8 kg
Capacity x Readability (kg)	/	3.2 x 0.0001	6.2 x 0.0001	8.2 x 0.001
Capacity x Readability (ct)	3100 x 0.5	16000 x 0.5	31000 x 0.5	41000 x 5
Capacity x Readability (oz)	21.870 x 0.005	112.880 x 0.005	218.700 x 0.005	289.25 x 0.05
Capacity x Readability (ozt)	19.930 x 0.005	102.880 x 0.005	199.335 x 0.005	263.60 x 0.05
Capacity x Readability (dwt)	398.7 x 0.1	2057.6 x 0.1	3986.7 x 0.1	5270 x 1
Capacity x Readability (lb)	1.3670 x 0.0005	7.0550 x 0.0005	13.6685 x 0.0005	18.080 x 0.005
Capacity x Readability (grn)	9570 x 2	49380 x 2	95680 x 2	126540 x 20
Stabilization Time (s)	1			
Construction	ABS plastic housing with 304 stainless steel (SST) pan			
Draftshield	No			
Calibration	SJX/E: User-selectable external span or linearity calibration SJX: Manual Internal Calibration / User-selectable external span or linearity calibration			
Tare Range	Full Capacity by subtraction			
Weighing Units**	g, kg, ct, oz, ozt, dwt, lb, grn, hkt, sgt, twt, tical, tola			
Application Modes**	Weighing, Parts Counting, Percent Weighing, Check Weighing, Totalization, Density Determination, Display Hold			
Typical Battery Life	120 hours	80 hours	120 hours	120 hours
Specified Temperature	10°C (50°F) to 40°C (104°F)			
Storage Conditions	-20°C (-4°F) to 55°C (131°F) at 10% to 90% relative humidity, non-condensing			
Communication	RS232, USB Host, USB Device, Ethernet or Bluetooth*** (all available as accessory)			
Display Type	Backlit LCD: 6-digit 7-segment with white LED backlight			
Display Size	20 mm digits / 0.78 in			
Pan Size (W x D)	170 x 140 mm / 6.7 x 5.5 in			
Scale Dimensions - SJX/E (W x D x H)	202 x 224 x 54 mm / 8 x 8.8 x 2.1 in			
Scale Dimensions - SJX (W x D x H)	202 x 230 x 68 mm			
Shipping Dimensions - SJX/E	300 x 250 x 86 mm / 11.8 x 9.8 x 3.4 in			
Shipping Dimensions - SJX (W x D x H)	300 x 250 x 129mm			
Net Weight (SJX/E)	1 kg / 2.2 lb			
Net Weight (SJX)	1.6 kg			
Shipping Weight (SJX/E)	1.5 kg / 3.3 lb			
Shipping Weight (SJX)	2.2 kg			

NTEP and Measurement Canada Approved models

Model	SJX323N/E	SJX622N/E	SJX1502N/E	SJX6201N/E
Capacity	320 ct (64 g)	620 g	1500 g	6200 g
Readability	0.01 ct (0.001 g) or 0.1 ct (0.01 g)	0.1 g	0.01 g or 0.1 g	1 g
Verification Interval e	0.1 ct (0.01 g)	0.1 g	0.1 g	1 g
Class	II	III	II	III
Span Calibration Mass*	50 g	300 g	1.5 kg	5 kg
Linearity Calibration Mass	30, 60 g	300, 600 g	1 kg, 1.5 kg	3 kg, 6 kg
Capacity x Readability (kg)	/	/	1.5 x 0.00001 or 1.5 x 0.0001	6.2 x 0.001
Capacity x Readability (ct)	320 x 0.01 or 320 x 0.1	3100 x 0.5	7500 x 0.5	31000 x 5
Capacity x Readability (oz)	2.2575 x 0.0005	21.870 x 0.005	52.910 x 0.005	218.70 x 0.05
Capacity x Readability (ozt)	2.0575 x 0.0005	19.935 x 0.005	48.225 x 0.005	199.35 x 0.05
Capacity x Readability (dwt)	41.153 x 0.001 or 41.15 x 0.01	398.7 x 0.1	964.52 x 0.01 or 964.5 x 0.1	3987 x 1
Capacity x Readability (lb)	/	1.3670 x 0.0005	3.3070 x 0.0005	13.670 x 0.005
Capacity x Readability (grn)	987.6 x 0.2	9568 x 2	23148 x 2	95680 x 20
Stabilization Time (s)	1.5	1	1.5	1
Construction	ABS housing & Stainless steel pan			
Draftshield	Yes	No		
Calibration	User-selectable external span or linearity calibration			
Tare Range	Full Capacity by subtraction			
Weighing Units**	g, kg, ct, oz, ozt, dwt, lb, grn			
Application Modes**	Weighing, Parts Counting, Percent Weighing, Check Weighing, Totalization, Density Determination, Display Hold			
Typical Battery Life	80 hours			
Specified Temperature Range	Class II: 10°C (50°F) to 30°C (86°F) Class III: 10°C (50°F) to 40°C (104°F)			
Storage Conditions	-20°C (-4°F) to 55°C (131°F) at 10% to 90% relative humidity, non-condensing			
Communication	RS232, USB Host, USB Device, Ethernet or Bluetooth*** (all available as accessory)			
Display Type	Backlit LCD: 6-digit 7-segment with white LED backlight			
Display Size	0.78 in / 20 mm digits			
Pan Size (W x D)	Ø93 mm / 3.7 in	Ø120 mm / 4.7 in	170 x 140 mm / 6.7 x 5.5 in	
Scale Dimensions (W x D x H)	202 x 222 x 103 mm / 8 x 8.7 x 4.1 in	202 x 224 x 54 mm / 8 x 8.8 x 2.1 in		
Shipping Dimensions (W x D x H)	300 x 250 x 129mm / 11.8 x 9.8 x 5.1 in	300 x 250 x 86 mm / 11.8 x 9.8 x 3.4 in		
Net Weight	2.2 lb / 1 kg			
Shipping Weight	3.3 lb / 1.5 kg			

COMPLIANCE

Compliance to the following standards is indicated by the corresponding mark on the product.

Mark	Standard
	This product complies with the EU Directives 2011/65/EU (RoHS), 2014/30/EU (EMC), 2014/35/EU (LVD) and 2014/31/EU (NAWI). The EU Declaration of Conformity is available online at www.ohaus.com/ce .
	This product complies with the EU Directive 2012/19/EU (WEEE) and 2006/66/EC (Batteries). Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. For disposal instructions in Europe, refer to www.ohaus.com/weee .
	AS/NZS 61000.6.1, AS/NZS 61000.6.3
	CAN/CSA-C22.2 No. 61010-1 UL Std. No. 61010-1

FCC Compliance Statement:

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
 - Increase the separation between the equipment and receiver.
 - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
 - Consult the dealer or an experienced radio/TV technician for help.
- Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

ISED Canada Compliance Statement:

This Class B digital apparatus complies with Canadian ICES-003.

INSTRUCTION MANUAL

For more detailed information about the balance, please refer to the product instruction manual. The instruction manual can be obtained either by scanning the QR code or by going to the Scout SJX product page to download it.

Go to www.ohaus.com to register your product's warranty.

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