Shaker Orbital SHHSMPDG Instruction Manual

EN

ES

FR

DE

IT

CZ

DA

NL

HU

PL

PT

SV

KR

CN

Table of Contents

1.	INTRODUCTION	2
1.1.	Safety Information	2
1.2.	Intended Use	2
1.3.	Package Contents	2
1.4.	Installation	3
1.5.	Overview	3
1.6.	Control Panel.....	4
2	OPERATION.....	5
2.1	Operation Instruction	5
2.2	Pulse time.....	5
3	MAINTENANCE.....	6
3.1	Troubleshooting.....	6
3.2	Service Information	6
4	TECHNICAL DATA.....	6
5	COMPLIANCE	7

1. INTRODUCTION

This manual contains installation, operation and maintenance instructions for the Ohaus shaker. Please read the manual completely before using.

1.1. Safety Information

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results.

WARNING For a hazardous situation with medium risk, possibly resulting in severe injuries or death if not avoided.

CAUTION For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or minor or medium injuries if not avoided.

ATTENTION For important information about the product. May lead to equipment damage if not avoided.

NOTE For useful information about the product.

Warning Symbols

General hazard

Electrical shock hazard

Safety Precautions

Please read the entire instruction manual before operating the unit.

WARNING: DO NOT use the UNIT in a hazardous atmosphere or with hazardous materials for which the unit was not designed. In addition, the user should be aware that the protection provided by the equipment might be impaired if used with accessories not provided or recommended by the manufacturer, or used in a manner not specified by the manufacturer.

Always operate unit on a level surface for best performance and maximum safety.

DO NOT lift the High-speed Microplate Shaker by the tray.

CAUTION! To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplug from the wall outlet. Disconnect unit from the power supply prior to maintenance and servicing.

Spills should be removed promptly

DO NOT immerse the unit for cleaning.

DO NOT operate the unit if it shows signs of electrical or mechanical damage.

Earth Ground – Protective Conductor Terminal.

Alternating Current

1.2. Intended Use

The Ohaus shaker is intended for general laboratory use. Safety cannot be guaranteed if used outside of the intended use.

1.3. Package Contents

- Microplate Shaker
- Power cord
- Instruction Manual

1.4. Installation


1. **Important:** choose the space on your bench where the Ohaus High Speed Microplate Shaker will stay. Because of the amount of energy dissipated, the unit is held onto your lab bench with 12 rubber feet (suction cups). The unit is VERY difficult to remove from work surface and should be done with caution. Use paper under the feet to keep them from attaching to the bench while determining the final location. When final location is determined, remove the paper.
2. Be sure the power switch is in the off position.
3. The Ohaus High Speed Microplate Shaker comes with a connector plug that is inserted into the IEC connector on the back of the unit first, and then it can be plugged into a properly grounded outlet. The 120v unit plugs into a 120 volt, 50/60 Hz source. The 230v unit plugs into a 230 volt, 50/60 Hz source.
4. Your Ohaus High Speed Microplate Shaker is ready for use.

1.5. Overview

Model:	SHHSMPDG
Controls:	See section 1.6
Overall dimensions (L x W x H):	15.5 x 12 x 12.8" (39.4 x 30.5 x 32.5cm)
Tray dimensions (L x W):	12.25 x 11" (31.1 x 27.9cm)
Power:	120V, 50/60Hz / 230V, 50/60Hz
Power consumption:	120W
Fuses:	5mm x 20mm, 5 amp quick acting (2 required)
Speed range:	600 to 2500rpm
Timer:	1 second to 160 hrs
Orbit:	.014" (3.6mm)
Net Weight:	48.1 lb / 21.8 kg
Ship weight:	50lbs / 22.7kg
Capacity:	7 lbs / 3.2 kg or 48 microplates in stacks up to 5" / 12.7 cm high
Communication:	RS-232 Interface

1.6. Control Panel

The front panel of the High Speed Microplate Shaker contains all the switches, controls and displays needed to operate the unit.

- A. **On/off rocker switch:** turns main power on/off.
- B. **Speed display:** Displays the speed of the shaker.
- C. Up/down arrows for set-point control.
- D. **Time display:** Displays accumulated time (continuous mode) or how much time is remaining (timed mode). The display range is from 0 to 9,999 minutes in one (1) second increments. The display will indicate minutes and seconds until the timer reaches 99 minutes and 59 seconds (99:59), then the display will automatically display minutes up to 9,999.
- E. Up/down arrows for set-point control.
- F. **Run/stop buttons:** Activates shaking.
- G. **Pulse on/off buttons:** Activates pulse mode.

2 OPERATION

2.1 Operation Instruction

For how to operate, please follow the following instructions.

- Flip up the four quick release latches which will release top tray from the 4 posts. DO NOT USE without both top and bottom tray foam in place.
- Your Digital High Speed Micro Plate Shaker is designed to hold up to 7 pounds, which will accommodate deep well blocks or stacked microplates up to 5" high. To avoid excessive noise and potential block damage, well blocks should always have space between them and should never touch each other.
- Sample containers should be placed on bottom tray foam inserted and centered on the tray to achieve optimal performance.
- To place top tray on unit, place the short side to the front engaging the 4 posts. When top tray is lined up, press down on the center of the tray and engage by pushing the quick release latches down into the locked position. Repeat for all corners. When top plate is locked in place, tighten the knobs on the four adjustment blocks. The High Speed Micro Plate Shaker is ready to operate. When changing samples, disengage the four quick release latches to remove the top tray. Do not loosen the knobs on the adjustment blocks. Loosen and move the four adjustment blocks only when the sample heights have changed.
- Turn the Power on by depressing the left side of the rocker switch on the bottom left side of the control panel. LED displays for Speed and Timer will light.
- Set speed by pressing the up/down arrows to the right of the speed display until you reach the desired speed. Adjust Time using the up/down arrows to the right of the timer display. Press the run button. The unit will run for the programmed time. Timer LED will display time remaining.
- To run in the un-timed mode, set the time to zero (0) using the up/down arrows to the right of the timer display. Press the run button and the unit will run until you press the stop button. Timer LED will display accumulated time.
- To run in pulse mode unit must not be running. Press the "ON" button above the word "PULSE". To stop the pulse feature press the "OFF" button above the word "PULSE". Pulse mode is factory programmed for 2 seconds on 1 second off. The pulse "ON" and pulse "OFF" times can be adjusted between 1 and 59 seconds in 1 second intervals.

2.2 Pulse time

Pulse "ON" time

- A. Turn power off
- B. Turn power on while holding pulse "ON" button. The LED display for "SPEED" shows current pulse on time.
- C. Change the pulse on time using the up/ down arrows to the right of the "SPEED" display.
- D. When the power is turned off, the value on the display is stored as the new pulse on time.

Pulse "OFF" time

- A. Turn power off
- B. Turn power on while holding pulse "OFF" button. The LED display for "TIMER" shows current pulse off time.
- C. Change the pulse off time using the up/down arrows to the right of the "TIMER" display.
- D. When the power is turned off, the value on the display is stored as the new pulse off time.

3 MAINTENANCE

The Ohaus High Speed Microplate Shaker is built for long, trouble-free, dependable service. No lubrication or other technical user maintenance is required. It needs no user maintenance beyond keeping the surfaces clean. The unit should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes. Spills should be removed promptly. Do not use a cleaning agent or solvent on the front panel which is abrasive or harmful to plastics, nor one which is flammable. Always ensure the power is disconnected from the unit prior to any cleaning. If the unit ever requires service, contact your Ohaus representative.

3.1 Troubleshooting

The following table lists common problems and possible causes and remedies. If the problem persists, contact OHAUS or your authorized dealer.

Error	Cause	Solution
Unit fails to power on	Missing or blown fuse	Add or replace fuse as necessary. If problem persists, please contact your OHAUS representative for repair.
Unit is excessively noisy	Rattling or ticking sounds may indicate a loose screw in the tray	Ensure that tray is secured tightly. If problem persists, please contact your OHAUS representative for repair.
Unit vibrates abnormally	Suction cup feet	Firmly press down on the four (4) corners of the unit creating a strong suction to the work surface (DO NOT place on bench mat). If problem persists, please contact your OHAUS representative for repair.

3.2 Service Information

If the troubleshooting section does not resolve or describe your problem, contact your authorized OHAUS service agent. For service assistance or technical support in the United States call toll-free 1-800-672-7722 ext. 7852 between 8:00 AM and 5:00 PM EST. An OHAUS product service specialist will be available to provide assistance. Outside the USA, please visit our website, www.ohaus.com to locate the Ohaus office nearest you.

Serial Number: _____

Date of Purchase: _____

Supplier: _____

4 TECHNICAL DATA

Operating Conditions: Indoor use only.

Temperature: 4 to 40°C (39.2 to 104°F)

Humidity: 20% to 80% relative humidity, non-condensing

Altitude: 0 to 6,562 ft (2000 M) above sea level

Non-Operating Storage:

Temperature: -20 to 65°C (-4 to 149°F)

Humidity: 20% to 80% relative humidity, non-condensing

Installation Category II and Pollution Degree 2 in accordance with IEC 664

5 COMPLIANCE

Compliance to the following standards is indicated by the corresponding mark on the product.

Mark	Standard
	This product complies with the applicable harmonized standards of EU Directives 2011/65/EU (RoHS), 2014/30/EU (EMC) and 2014/35/EU (LVD). The EU Declaration of Conformity is available online at www.ohaus.com/ce .
	This product complies with the EU Directive 2012/19/EU (WEEE). Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. For disposal instructions in Europe, refer to www.ohaus.com/wEEE .
	EN 61326-1
	CAN/CSA-C22.2 No. 61010-1, CAN/CSA-C22.2 No. 61010-2-051 UL 61010-1, UL 61010-2-051

ISED Canada Compliance Statement

This Class A digital apparatus complies with Canadian ICES-003.

ISO 9001 Registration

The management system governing the production of this product is ISO 9001 certified.

FCC Supplier Declaration of Conformity

Unintentional Radiator per 47CFR Part B

Trade Name: OHAUS CORPORATION

Model or Family identification: SHHSMPDG

Party issuing Supplier's Declaration of Conformity:

Ohaus Instruments (Changzhou) Co., Ltd.

2F, 22 Block, 538 West Hehai Road, Xinbei District, Changzhou

Jiangsu 213022,

China

Phone: +86 519 85287270

Responsible Party – U.S. Contact Information:

Ohaus Corporation

7 Campus Drive, Suite 310

Parsippany, NJ 07054

United States

Phone: +1 973 377 9000

Web: www.ohaus.com

FCC Compliance Statement:

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

LIMITED WARRANTY

OHAUS products are warranted against defects in materials and workmanship from the date of delivery through the duration of the warranty period. During the warranty period OHAUS will repair, or, at its option, replace any component(s) that proves to be defective at no charge, provided that the product is returned, freight prepaid, to OHAUS.

This warranty does not apply if the product has been damaged by accident or misuse, exposed to radioactive or corrosive materials, has foreign material penetrating to the inside of the product, or as a result of service or modification by other than OHAUS. In lieu of a properly returned warranty registration card, the warranty period shall begin on the date of shipment to the authorized dealer. No other express or implied warranty is given by OHAUS Corporation. OHAUS Corporation shall not be liable for any consequential damages.

As warranty legislation differs from state to state and country to country, please contact OHAUS or your local OHAUS dealer for further details.

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SAP P/N 30570467A P/N 780166-00Rev2 © 2020 Ohaus Corporation, all rights reserved