

ULTRASONIC WATER BATHS

(AMBIENT +5°C TO +50°C)

Image: UB-420, UB-410 and UB-405



Powersonic Digital Ultrasonic Water Baths

- Adjustable temperature range up to 50°C
 - Accurate digital PID temperature controller
 - Durable ABS moulded plastic body
 - Seamless stainless steel tank with radial corners
 - Adjustable digital timer from 0 to 99 minutes
 - Stainless steel baskets and moulded lid supplied standard
 - Built in drain valve
-

PRODUCT DETAILS

Three powerful and reliable models available in the digital Powersonic range. Each model has a built in microprocessor, built in thermometer, 0 to 99 minute timer with alarm to indicate cycle end, built in drain valve and perforated stainless steel baskets.

The units are made with moulded plastic which makes them durable and long lasting.

SPECIFICATIONS

Model	UB-405	UB-410	UB-420
Capacity	5 litres	10 litres	20 litres
Bath Dimensions (LxWxD mm)	300x155x150	300x240x150	500x300x150
External Dimensions (WxDxH mm)	500x285x255	510x275x275	720x450x275
Temperature Range	Ambient to 50°C		
Temperature Control	Digital PID		
Transducer	Industrial BLT type		
Frequency	40kHz		
Housing Material	ABS moulded plastic		
Power (W)	350	500	700
Weight (kg)	7	9	14

Thermoline Scientific have been manufacturing and distributing high quality laboratory and scientific testing equipment since 1970. Over this time, Thermoline has grown to be a leading brand in the science industry, with our equipment being used in small and large Hospitals, Universities and Research Laboratories across Australia and the Asia Pacific region.

We're proud to say that we are 100% Australian owned and operated.

Head Office Phone (02) 9604 3911

Head Office Fax (02) 9725 1706

Email info@thermoline.com.au

Web www.thermoline.com.au