



Benchtop pH/Water Quality Analyzer

LAQUA

F-70/DS-70 Series

Intuitive and very easy-to-use touch panel operation

- pH**
- mV(ORP)**
- ION**
- Conductivity**
- Resistivity**
- Salinity**
- TDS**



Simply slide your finger across the screen to switch displays



Two channels can be displayed simultaneously



Color LCD touch panel display


| | | | |
|---|--|--|---|
| <p>NAVI 2CH USB PC PRT ID USP/EP/JP</p> <p>F-74 CH.1 pH ORP ION CH.2 COND RESI SAL TDS</p> | <p>NAVI 2CH USB PC PRT ID USP/EP/JP</p> <p>F-73 CH.1 pH ORP ION CH.2 pH ORP ION</p> | <p>NAVI 2CH USB PC PRT ID USP/EP/JP</p> <p>F-72 CH.1 pH ORP ION</p> | <p>NAVI USB PC PRT ID USP/EP/JP</p> <p>DS-72 COND RESI SAL TDS</p> |
|---|--|--|---|

*Accessories included : Electrode stand/Instruction manual/Quick manual/AC adapter/Cover (F-72/F-73/F-74/DS-72 only)

Full support for on-screen setting confirmation, maintenance information and troubleshooting tips guide you through trouble-free operation


Inspection Navigation

Easy navigation for main unit and electrode inspections. Various industrial standards (JIS, USP, EP, JP, CP) are also supported.




Troubleshooting Navigation

Reliable on-screen support if a problem occurs during calibration or measurement. The software has a user guide to resolve any operation difficulties.



Application Functions

Various industry standard methods are supported from measurement to result output. Conductivity measurements for pharmaceutical pure water guidelines of various countries are also supported.



Free Arm Electrode Stand

The free arm of the stand-alone electrode stand can be positioned in any direction, vertically or horizontally. The long-type electrode stand* with a telescopic stand is also provided for measurements with large beakers.

*Optional



Full-Range Functions for Validation and Usability

- Periodic inspection mode: JIS/Pharmacopeias/Digital Simulator (F-72/F-73/F-74)
- Full support for pharmaceutical pure water guidelines of various countries. (USP/EP/JP/CP) (F-74/DS-72)
- Customizable auto-hold function for calibration and measurement (F-72/F-73/F-74/DS-72)
- Simultaneous connection to a GLP/GMP compatible printer and PC
- Digital memory: Up to 2,000 sets of measurement data can be recorded (F-71/F-74BW/DS-71:999)
- USB-PC communication *(all models) and USB memory (F-72/F-73/F-74/DS-72)
- Multi-language support (Japanese, English, Chinese, Korean) (F-72/F-73/F-74/DS-72)
- FDA21CFR Part 11 (please ask for a quotation)

Custom LCD display

