



# Compact Scale

*Portable Electronic Balances*



Ohaus Compact Scales are ideal when cost matters and dependable results are a must. These uniquely designed scales have a wide variety of applications, more than any other scale in their class. The three portable models range in capacities from 200 g to 5000 g, and are perfect for use in light production, in the field, the kitchen, the office, and more. Easy operation, long battery life, and great performance make Compact Scales a truly remarkable value.

- Sleek, low profile design
- Large high contrast display
- Energy-saving auto shut-off feature
- Battery or AC powered
- Rugged plastic construction

# CS Series Portable Electronic Balances

## Specifications

Capacity (g)	200	2000	5000
Readability (g)	0.1	1	2
Repeatability (Std. Dev.) (g)	0.1	1	2
Linearity (g)	0.1	1	2
Weighing Units	g, lb:oz		
Tare Range	to capacity by subtraction		
Stabilization Time	3 seconds		
Optimum Operating Temp.	18° - 25°C		
Power	3 AA batteries (included) AC adapter (optional accessory)		
Auto shut-off	5 minutes no activity (battery operation only)		
Calibration	digital calibration from keypad		
Display (cm)	LCD (11)		
Pan Size (cm)	146 x 133		
Dimensions LxWxH (cm)	193 x 135 x 39		
Net Weight (kg)	0.47		
Item Number	CS200	CS2000	CS5000

## Standard Features


Low battery indicator, overload and underload indicators, auto shut-off, rugged plastic construction, tactile feedback keypad.

## Accessories

1050 ml capacity bowl, AC adapter, stainless steel pan cover, calibration masses.

## Simple To Order

Contact your dealer for product availability as well as warranty information.



**Ohaus UK Ltd.**  
64 Boston Road  
Beaumont Leys  
Leicester LE4 1AW  
ENGLAND  
Tel.: +44 116 234 5075  
Fax: +44 116 234 6706  
e-mail: [uksales@ohaus.com](mailto:uksales@ohaus.com)

**Ohaus Europe & MEA**  
Im Langacher  
CH-8606 Greifensee  
Switzerland  
Tel.: +41-1-944 33 66  
Fax: +41-1-944 32 30  
e-mail: [EuropeSales@Ohaus.com](mailto:EuropeSales@Ohaus.com)

*With offices in Mexico, Switzerland,  
England, Germany, France, Spain, Italy,  
Poland, Russia, Japan, Korea and China*

\* ISO 9001  
Registered Quality Management System

[www.ohaus.com](http://www.ohaus.com)

80774070 © OHAUS Europe GmbH

