

## WRESTLING REFRACTOMETER

Refractometers have been proven time and time again to be more accurate than hydrometers or reagent test strips for testing urine specific gravity. The Palm Abbe is fast, convenient, and easy to use. Simply place a drop or two of urine in the stainless-steel sample well, press a button, and instantly read urine specific gravity on the large LCD display.

Model: PA201X Part No.: 093

Scale #	Unit of Measure	Range	Resolution	Precision (+/-)
093	Urine Specific Gravity (D20/20)	1.0000 to 1.0450	0.0001	0.0005

	MISCO Palm Abbe	Competitor
Prism Material	SAPPHIRE (Comparable to Diamond)	GLASS
Detector Array Elements	1,024 Elements	128 Elements
Detector Resolution	3,256 ppi	400 ppi
Stainless Steel Well	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Number of Scales	Up to Five	One
Display	24 Alphanumeric Characters x2 Lines	Four Digits
Display Backlight	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Evaporation Cover	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Optional Rubber Armor	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Made In	U.S.A	Asia

Made in the USA! - Copyright © 2013 MISCO