

PREMIUM BLANKET WARMERS

(AMBIENT +5°C TO +60°C)

Images: TBW-520F-1-GD and TBW-1100F-2-GD



These premium Blanket Warmers models have been designed to accommodate medical centres and hospitals where blankets are required to be stored safely in a warm and clean environment.

- Stainless steel internal liner which is easy to clean and is corrosive resistant
- Digital temperature display with high temperature alarm
- Gravity convection or internal forced air flow (210 litre model)
- 13mm cable/monitoring port hole complete with rubber plug
- Adjustable shelves
- Glass door/s for easy viewing
- Internal light on selected models
- Castors fitted for easy of mobility on selected models (excludes 210 litre model)



PRODUCT DETAILS

Thermoline has four key blanket warmer models all with precise temperature control between ambient +5°C to +60°C.

For added safety, the digital temperature control incorporates an automatic over temperature alarm which in the case of an over temperature condition, the alarm will be activated and will shut down power to the heating elements. This occurs as well as emitting a visual signal to warn the operator of an alarm condition. Connection to a building management system (BMS) is available as standard on all models.

Thermoline's blanket warmers have internal fan forced air circulation, energy efficient LED lights plus easy view glass self-closing and lockable doors.

All models are supplied complete with heavy duty castors to allow cabinets to be easily moved when the area around cabinet requires cleaning.

A 13mm Ø cable/monitoring port hole with removable rubber plug is available on the side of each model. This allows easy access for a remote temperature measurement device.

Safety Features

Deviation high alarm on temperature controller (Failsafe)



Image: TBW-210F-1-GD

SPECIFICATIONS

Model	TBW-210F-1-GD	TBW-360F-1-GD	TBW-520F-1-GD	TBW-1100F-2-GD
Internal Dims. (WxDxH mm)	500x505x735	590x550x990	590x550x1460	1320x550x1460
External Dims. (WxDxH mm)	625x595x985	740x730x1560	740x730x2010	1470x730x2010
Temp. Range	Ambient +5°C to 60°C			
Capacity (litres)	210	360	520	1100
Door Type	Glass Door/s			
Shelves	3 Shelves		4 Shelves	4 Levels
Electrical (W / V)	415 / 240	450 / 240	450 / 240	900 / 240
Weight (kg)	60	120	180	380

Head Office Phone (02) 9604 3911

Head Office Fax (02) 9725 1706

Email hello@thermoline.com.au

Web www.thermoline.com.au



Thermoline Scientific have been manufacturing and distributing high quality laboratory and scientific testing equipment since 1970. Over this time, Thermoline has grown to be a leading brand in the science industry, with our equipment being used in small and large Hospitals, Universities and Research Laboratories across Australia and the Asia Pacific region.

We're proud to say that we are 100% Australian owned and operated.