

Oxford BTD

Bottle Top Dispenser

Made of highly chemical-resistant materials, such as PTFE and PFA, the Oxford Lab Products Bottle-Top Dispenser is the ideal choice for dispensing bases, acids and solvents.

Eco-friendly - while discharging the tubes in the re-circulation mode, in fast priming, the reagent is being redirected back into the bottle, eliminating any waste of precious reagents, ensuring safety and savings.

Unique piston mechanism - cleaning of piston and cylinder without affecting the calibration.

Dual graduation marks - readability of the volume from both sides
Springless valve design - leads to smooth functioning and high chemical resistance.

Easy handling - discharge tube with 360° rotation, makes the bottle label visible all the time / telescoping filling tube compatible with different sizes of reagent bottles / ergonomic grip.

Autoclavable at 120 degrees Celsius for 20 minutes.

Easy to calibrate - complies with ISO 8655-5 standard.

Adapters- comes with 4 bottle adapters for the most common bottle sizes.



Model#	Description	Volume
AD-2.5	AccuTop Bottle Top Dispenser, Without Re-circulation	0.25-2.5 ml
AD-05	AccuTop Bottle Top Dispenser, Without Re-circulation	0.5-5 ml
AD-10	AccuTop Bottle Top Dispenser, Without Re-circulation	1-10 ml
AD-30	AccuTop Bottle Top Dispenser, Without Re-circulation	2.5-30 ml
AD-60	AccuTop Bottle Top Dispenser, Without Re-circulation	5-60 ml
ADR-2.5	AccuTop Bottle Top Dispenser, With Re-circulation	0.25-2.5 ml
ADR-05	AccuTop Bottle Top Dispenser, With Re-circulation	0.5-5 ml
ADR-10	AccuTop Bottle Top Dispenser, With Re-circulation	1-10 ml
ADR-30	AccuTop Bottle Top Dispenser, With Re-circulation	2.5-30 ml
ADR-60	AccuTop Bottle Top Dispenser, With Re-circulation	5-60 ml

10 Year Warranty

For more details, visit us at www.oxfordlp.com