

HVL-CS

Digital Chair Scale



- **Accurate Patient Weighing**
- **Comfortable easy to clean seat**
- **Easy to read display**
- **Convenient push handles**

AND
A&D Australasia
www.andweighing.com.au





Features

Ergonomically designed

Minimal effort required by patient and staff

Simple operation

2 button keypad

Easy to clean chair shell

Ensures maximum hygiene

Long battery life

600 hrs for standard alkaline cells

Auto power function

Scale switches off after 3 minutes of inactivity preserving battery life

Highly portable

Simple to move using the frictionless wheels

Locking castors

To ensure patient safety when weighing

High visibility LCD

Easy to view from any angle

Push handles

Convenient for staff when moving the scale

Australian made

Quality manufacturing (ISO-9001)



Specifications

Capacity	150kg
Minimum Division	20g
Weight Display	Liquid Crystal
Power	6 x 'D' Cell batteries or AC adaptor (optional)
Chair Shell	White high strength moulded plastic
Base	Grey die cast aluminium spider
Wheels	5 x 125 mm diameter
Safety	All wheels fitted with locking brakes
Loadcell	Allows 200% overload without failure

Dimensions

	(W) x (D) x (H) (mm)
Overall Dimensions	750 x 630 x 880
Packed Dimensions	930 x 670 x 950

Specifications are subject to change for improvement without notice.
Photographs shown are for illustration purposes only.