





PAL-BX/RI features a dual scale. Brix (wide range 0.0-93.0%) and Refractive Index (resolution 0.0001).

Digital Hand-held "Pocket" **Dual Scale Refractometer** 

Cat.No.3852

(C) ATAGO®

### Measurement Method



1. Place the sample onto the prism surface.



2. Press the START key.



3. The measurement value will be displayed in 3 seconds.

4. Press the START key for 2 seconds while the Brix (%) is indicated, to display the Refractive Index and temperature.

Press the START key (one time) to return to the Brix (%) scale.

### Example of display

Display of the measurement value "1.3425" and "20°C".







The measurement value for Refractive Index will scroll from right to left.













Actual size

# **Specifications**

## **Features**

- Requires only a few drops of sample
- Measurement in 3 seconds
- Ergonomically designed for easy, one-handed operation
- Light & Compact, 100g
- Designed with a strap holder option
- Highly water resistant (IP65)
- Convenient storage case
- Battery power indicator

Brix : 0.0 to 93.0% (9.0 to 99.9°C) Refractive Index (RI) : 1.3306 to 1.5284 (5.0 to 45.0°C)
Brix: 0.1% Refractive Index: 0.0001 Temperature: 0.1°C
Brix: ±0.1% Refractive Index: ±0.0003 (Water at 20°C) Temperature: ±1°C
10 to 100°C (Brix only)
10 to 40°C
AAA alkaline battery × 2
Approx. 11,000 measurements (using new alkaline batteries)
IP65 Water resistant
55(W) × 31(D) × 109(H)mm, 100g (main unit only)

<sup>\*</sup> The displayed measurement value for Refractive Index is not temperature compensated. The reading is the Refractive Index of the sample at the temperature displayed.

All ATAGO refractometers are designed and manufactured in Japan.



ATAGO products comply with HACCP, GMP, and GLP system standards.



Headquarters: The Front Tower Shiba Koen, 23rd Floor 2-6-3 Shiba-koen, Minato-ku, Tokyo 105-0011, Japan TEL: 81-3-3431-1943 FAX: 81-3-3431-1945 overseas@atago.net http://www.atago.net/



TEL: 1-425-637-2107 customerservice@atago-usa.com

TEL: 91-22-2833-8038 / 8076 customerservice @ atago-india.com

ATAGO (THAILAND) Co., Ltd.
TEL: 662-982-8718-9 customerservice@atago-thailand.com

CATAGO BRASIL Ltda.
TEL:55 16 3916-6000 customerservice@atago-brasil.com

ATAGO ITALIA s.r.I.
TEL: 39 2 36557267 customerservice@atago-italia.com

ATAGO CHINA Guangzhou Co., Ltd.
TEL: 86 20-38106065 info@atago-china.com



\* Specifications and appearance are subject to change without notice.