

PREMIUM FLUID WARMERS

(AMBIENT +5°C TO +40°C)

Image: TFW-520-2-GD



Thermoline Premium Fluid Warmers

- Internal forced air flow
- Digital temperature display with high and low temperature alarms
- Stainless steel internal liner which is easy to clean and is corrosive resistant
- Internal light
- Cable/monitoring port hole
- Double glazed glass doors
- Self closing doors with key locks
- Adjustable PVC coated removable shelves
- Castors fitted for easy moving

Head Office Phone (02) 9604 3911
Head Office Fax (02) 9725 1706

Email info@thermoline.com.au
Web www.thermoline.com.au

 **Thermoline**
S C I E N T I F I C
H I G H E R S T A N D A R D S I N S C I E N C E

PRODUCT DETAILS

Two premium fluid warmer models are available, a 520 litre model with a single front door and a 1100 litre model with double front opening doors. Both are designed to accommodate the requirements of medical centres and hospitals where fluids are required to be stored safely.

Temperature is precisely controlled electronically between ambient +5°C to +40°C. For added safety the digital temperature control incorporates an automatic over temperature alarm which in the case of a over temperature condition the alarm will be activated and shut down power to the heating elements as well as emit a visual signal to warn the operator of an alarm condition. Connection to a building management system (BMS) is available as an optional extra.

Cabinet construction is grey Colorbond exterior with interior finish of 304 stainless steel. Insulation is high efficiency injected polyurethane foam. Both models are supplied standard with double glazed glass self closing doors with key locks on the doors for security.

Optional Accessories

Additional PVC coated shelves

Remote alarm contacts

Safety Features

Deviation high alarm on temperature controller (Failsafe)

Manual resetting element safety

SPECIFICATIONS

Model	Internal Dimensions (WxDxH mm)	Capacity (litres)	External Dimensions (WxDxH mm)	Shelves	Electrical	Weight (kg)
TFW-520-1-GD	580x560x1460	520	720x720x2040	4 levels	3A	180
TFW-1100-2-GD	1320x560x1460	1100	1470x770x2040	4 levels	5A	380



Thermoline Scientific have been manufacturing and distributing high quality laboratory and scientific testing equipment since 1970. Over this time, Thermoline has grown to be a leading brand in the science industry, with our equipment being used in small and large Hospitals, Universities and Research Laboratories across Australia and the Asia Pacific region.

We're proud to say that we are 100% Australian owned and operated.